

NEW!

SIMPLY SMARTER CRYOSURGERY



INTRODUCING



**CRYOSURGERY WITHOUT
THE GUESSWORK**

- ❄️ **Nitrous Oxide OR Carbon Dioxide Gas**—First and only cryo product that offers the flexibility to use either cryogen based on your patients' needs.
- ❄️ **Advanced Treatment Wand**—Now you can treat your patients in any direction, at any angle, on any part of the body.
- ❄️ **Preset Dispense Times**—Choose between 5-second or 10-second treatment times and minimize the risk of under- or over-treating.
- ❄️ **Gas Level Indicator**—Our new lighted indicator ensures you will know when to order a new tank.



CRYOSURGERY JUST GOT A WHOLE LOT SMARTER!


For the first time ever, you can choose between two cryogens—Nitrous Oxide (-89°C) or Carbon Dioxide (-79°C)—to get great first-time outcomes and minimize damage to surrounding skin.

N₂O is an ideal cryogen for most common lesions such as verruca vulgaris, angiomas, seborrheic keratosis and more.*

What cryogen is best for cosmetic lesions such as skin tags, age and sun spots on the hands and face? The answer is CO₂ which provides a more superficial treatment and minimizes your patients' downtime and recovery.



FREEZING POINT COMPARISON



Item No.	Product Code	CryoLab® Medical Product Description
Please contact CryoConcepts, LP for details 200-1000 MAR-078 (REV 3)	CryoLab Medical Device Unit	Kit Contains: 1 Canister of cryogen 1 N ₂ O OR CO ₂ 1 12V/20W Power Supply 1 90" Power Cord 1 Quick Start Guide 1 Instructions for Use One-year Limited Warranty
000-2730	Replacement N ₂ O Cylinder	Contains: 1 Canister of N ₂ O cryogen 1 Bleed Plug 1 Return Instructions 1 Return Label
000-2731	Replacement CO ₂ Cylinder	Contains: 1 Canister of CO ₂ cryogen 1 Bleed Plug 1 Return Instructions 1 Return Label



Contact your Prime Vendor TODAY!



Providing the best options the industry has to offer.

Jonathan Zarnowski
 Direct Line: 248-781-8700
 Email: JZarnowski@ELPROVIDERS.COM
 Website: ELPROVIDERS.COM