E-MY INTERIOR OFFICE THEFAULT BETTER OUTCOMES

EARLY DETECTION. RAPID DIAGNOSIS. B

LY DETECTION, RAPIS BIAGINOSIS, BETTER BUTCOMES.

# RLY DETECTION. RAPID DIAGNOSIS. BE



## **PRODUCT CATALOG**

EARLY DETECTION. RAPID DIAGNOSIS. BETTER OUTCOMES:

Y DETECTION PAPID DIAGNOSIS RETTER OUTCOMES EARLY DETECTION, RAPID BIAGNOSIS

ER OUTCOMES. EARLY DETECTION. RAPID DIAGNOSIS. BETTER

DIAGNOSIS: BETTER OUTCOMES.

RLY DETECTION. RAPID DIAGNOSIS. BETTER OUTCO

DIAGNOSIS. BETTER OUTCOMES.

# CTION. RAPID DIAGNOSIS. BETTER OUT

EARLY DETECTION, RAPID DIAGNOSIS, BETTER OUTCOMES.

EARLY DETECTION. RA

TER OUTCOMES. FARLY DETECTION. RAPID DIAGNOSIS. BETTER

TION RAPID DIAGNOSIS BETTER OUTCOMES

RLY DETECTION. RAPID DIAGNOSIS. BETTER OUTCO

DIAGNOSIS. BETTER OUTCOMES

CTION. RAPID DIAGNOSIS. BETTER OUT

RAPID DIAGNOSIS • EARLY DETECTION • BETTER OUTCOMES

# CTION. RAPID DIAGNOSIS. BETTER OUT

EARLY DETECTION PAPID DIAGNOSIS BETTER OUTCOMES

EARLY DETECTION. RAPID DIAGNOSIS. BETTER OUTCOMES.

BETTER OUTCOMES. EARLY DETECTION. RAPID DIAGNOSIS. BET

EARLY DETECTION. RAPID DIAGNOSIS. BETTER OUTCOMES

EARLY DETECTION DARID DIACNOSIS DET

LifeSign is a medical diagnostics company delivering high-quality rapid point-of-care testing products to the healthcare industry.

LifeSigns's extensive line of innovative products is used by clinicians worldwide to aid in the detection of medical conditions which include Infectious Diseases, Women's Health, Drugs of Abuse, and Cardiac Markers.

For pricing and purchasing, please reach out and contact:

**2**48-781-8700

**1-877-745-9943** 

www.elproviders.com

## **CARDIAC CARE**



LifeSign MI

For the Rapid Qualitative Detection of Cardiac Markers to Aid in the Diagnosis of Acute Myocardiac Infarction in Human Whole Blood, Serum or Plasma.

- Qualitative CK-MB, Myoglobin, Troponin I
- Enables prompt Emergency Room decisions
- Easy to use procedure
- Results in 15 minutes
- Reduces unnecessary testing and admissions

ORDERING INFORMATION		
Product	# OF TESTS	ITEM #
LifeSign MI Troponin I	20	60101
LifeSign MI Myoglobin/CK-MB/Troponin I	20	60201
LifeSign MI Cardiac Control Set	2 X 1 ml	LSC-301

## **UROLOGY**

#### Uricult®

The Preferred In-Office Urine Culture System.

- Easy to perform; Dip Incubate Read Results
- Results within 18-24 hours (Incubation)
- Positive cultures can be sent to Reference Lab for additional studies if needed
- Trademark of Orion Diagnostica

Product	# OF TESTS	ITEM #
Uricult CLED/EMB	10	1000
Uricult CLED/Polymixin/EMB	10	1002
Uricult CLED/Polymixin/MacConkey	10	1003





#### **Uricult® Control Kit**

Quality Control for Uricult Urine Culture Paddles. Each kit contains a positive and a negative control plus diluent.

ORDERING INFORMATION		
Product	# OF TESTS	ITEM #
Uricult Control Kit	2 x 99ml	200030

## **UROLOGY**



#### Status BTA

A Rapid Test for the Qualitative Detection of Bladder Tumor Associated Antigen in Human Urine.

- CLIA Waived
- Indicated for use as an aid in monitoring bladder cancer patients in conjunction with cystoscopy
- Accurate results with superior sensitivity & specificity
- Results in 5 minutes

ORDERING INFORMATION		
Product	# OF TESTS	ITEM#
Status BTA	25	39025
Status BTA Control Set	2 X 2 ml	200019

## **INFECTIOUS DISEASE**

#### Status Flu A and B

For the Qualitative Detection of Influenza Type A and Type B Antigens Directly from Nasal Swab, Nasopharyngeal Swab and Nasal Aspirate/Wash Specimens.

- CLIA Waived for Swab Specimens
- Flexible and flocked nasal swab for improved patient comfort and superior specimen collection
- Innovative flip test design with on board sample extraction
- Positive and Negative control swabs included
- Differentiates between Influenza A and B
- Individual pre-measured extraction reagent capsule for increased accuracy and ease of use

ORDERING INFORMATION		
Product	# OF TESTS	ITEM #
Status Flu A&B test kit	25	36025
Nasal Aspirate Kit	22	510AS
Status Flu A&B Control swabs	5 sets (neg+pos)	36022CS
Status Flu A&B nasal swabs refill	52	36052





#### Status Mono

For the Qualitative Detection of Infectious Mononucleosis Heterophile Antibodies in Whole Blood, Serum or Plasma.

- · FDA cleared with no age restrictions
- · CLIA Waived for Whole Blood
- Results within 8 minutes
- Requires only 25ul of whole blood or 10ul of serum/ plasma
- Superior 99% sensitivity specificity of 98.8%
- External positive and negative control included (item #68364 and 84M30)
- Safe and accurate plastic sample pipettes

ORDERING INFORMATION		
Product	# OF TESTS	ITEM #
Status Mono (WB/Serum/Plasma)	25	68364
Status Mono (WB/Serum/Plasma)	30	84M30
Status Mono (WB only)	30	84W30
Status Mono (WB only)	10	84W10
Status First Mono Control Set	2 X 1 ml	200014

## **INFECTIOUS DISEASE**

#### **Status Strep A Products**

For the Rapid Detection of Group A Streptococcal Antigen Directly from Throat Swab Specimens.

#### Status Strep A Flip

- CLIA Waived
- Patented single use reagent capsules
- Positive results in 5 minutes
- Sensitivity 96%, Specificity 99%
- No counting reagent drops-less steps for greater accuracy





## **Status Strep A Strip**

- CLIA Waived
- Sensitivity of 96.2%, specificity of 98.7%
- Results in 5-10 minutes
- · External positive and negative controls included

## Status Strep A Plus

- CLIA Waived
- Sensitivity of 96.2%, specificity of 98.7%
- Results in 5-10 minutes
- · External positive and negative controls included
- Two convenient resealable canisters holding 25 test strips each





## Status Accustrep Strep A

- Sensitivity of 96.2%, specificity of 98.7%
- Results in 5-10 minutes

ORDERING INFORMATION		
Product	# OF TESTS	ITEM #
Status Strep A Flip (cassettes)	25	34125
Status Strep A Strip (pouch)	30	34130
Status Strep A Plus (cannister-strips)	50	34250
Status Accustrep Strep A (cassettes)	25	34025

## **INFECTIOUS DISEASE**

## **Uni-Gold Recombigen HIV**

For Qualitative Detection of HIV-1/2 in Serum, Plasma and Whole Blood.

- · CLIA Waived for whole blood
- · Results in 10-12 minutes
- Sensitivity of 100%, Specificity of 99.7%
- Trademark of Trinity Biotech

ORDERING INFORMATION		
Product	# OF TESTS	ITEM #
Uni-Gold HIV	20	1206506
Uni-Gold HIV Control Set	2 X 0.5 ml	1206530



## **GI DISORDERS**



### Status H. pylori

For the Qualitative Detection of Anti-Helicobacter pylori Antibody in Human Whole Blood, Serum and Plasma.

- CLIA Waived for whole blood
- Requires only 25ul of whole blood or 10ul of serum or plasma
- Results in 10 minutes
- Sensitivity of 95.9% with overall agreement to biopsy histology

ORDERING INFORMATION		
Product	# OF TESTS	ITEM #
Status H. pylori	30	37030
Status H. pylori Control Set	2 X 1 ml	200013

#### Status IFOBT

Rapid Immunoassay for the Qualitative Detection of Immunochemical Fecal Occult Blood. (Human Hemoglobin in Human Fecal Specimens)

- CLIA Waived
- No dietary or medication restrictions
- · Increased patient compliance
- Results in 5 minutes
- Contains 30 test devices, 50 collection tubes, 50 mailing envelopes

#### www.statusifobt.com

ORDERING INFORMATION		
Product	# OF TESTS	ITEM #
Status IFOBT (complete kit)	30/50/50	38030-DTM
Status IFOBT Sample Collection Tube Only	50	38050-T
Status IFOBT Mailing Envelope Only	50	38050-M
Status IFOBT Mailers & Collection Tubes	25	38025-TM
Status IFOBT Device & Tubes	50	38050-DT
Status IFOBT Control Set	2 X 1 ml	200020



## **WOMEN'S HEALTH**

#### Status hCG

For the Detection of hCG in Serum or Urine.

- **CLIA Waived** for urine
- Serum and urine 25mIU/ml sensitivity
- Results in 3-5 minutes
- Minimal susceptibility to the Hook effect by hCG  $\beta$ core fragment



**Status hCG Strip**For the Detection of hCG in Urine.

- **CLIA Waived**
- Simple "dip and read" procedure 25mlU/ml sensitivity
- Results in 3-5 minutes
- Minimal susceptibility to the Hook effect by hCG  $\beta$ core fragment



ORDERING INFORMATION		
Product	# OF TESTS	ITEM #
Status hCG Urine/Serum	35	35035
Status hCG Urine	35	35135
Status hCG Strip Urine	35	35235
Status hCG Control Set	2 X 2.5ml	135-100A

## **DRUGS OF ABUSE**

One Step Immunoassay for the Qualitative Detection of Drugs of Abuse and their Metabolites in Urine.

#### Status DS

- Simple one-step procedure
- Accurate results with 97-99% correlation to GC/MS
- Results within 5-10 minutes
- Flexible combination of drugs up to a 10 panel cassette

ORDERING INFORMATION		
Product	# OF TESTS	ITEM #
Status DS TCA	35	20935
Status DS 8 Panel THC/OPI/COC/AMP/BZO/ BAR/TCA/PCP	25	28025
Status DS 10 Panel MET/THC/OPI/COC/AMP/ BZO/BAR/TCA/PCP/MTD	10	21010
Status DS 10 Panel MET/THC/OPI/COC/AMP/ BZO/BAR/TCA/PCP/MTD	25	21025
Status DS Control Set	2 X 5 ml	21000



## Status Toxi-Cup 13+4

- CLIA Waived
- Simple Procedure with easy to read results
- Results in 5 minutes
- Adulterants Included
- Temperature Indicator

ORDERING INFORMATION		
Product	# OF TESTS	ITEM #
Status Toxi-Cup 13+4: AMP/BAR/BUP/BZD/ COC/MAD/MDMA/MET/OPI/OXY/PCP/TCA/THC Adulterants - CRE/pH/SG/NBP	25	65C13+4
Status DS Control Set	2x5 ml	21000



# NICOTINE

#### **Status Nicotine**

Rapid Immunoassay for the Detection of Cotinine, Nicotine Metabolite in Urine.

- Simple one step procedure
- Results in 10 minutes
- 500ng/ml cut-off

ORDERING INFORMATION		
Product	# OF TESTS	ITEM #
Status DS Nicotine	35	21735
Status DS Control Set	2 x 5 ml	21000



## **ALCOHOL**

#### Q.E.D. A150®

For the Rapid Quantitative Determination of Alcohol in Saliva.

- CLIA Waived
- Simple two step procedure
- Accurate results; built-in QA spot<sup>™</sup> assures quality control

ORDERING INFORMATION		
Product	# OF TESTS	ITEM #
Q.E.D A150	10	31150



## **HYGIENE**



## Hygicult®

Rapid Hygiene Culture System.

- · Ready to use for on site hygiene monitoring
- Store at room temperature
- Hinged slide for convenient surface sampling
- AOAC approved
- Incubate at 36°C for 24 hours
- •Trademark of Orion Diagnostica

ORDERING INFORMATION		
Product	# OF TESTS	ITEM #
Hygicult™ TPC	10	68010
Hygicult™ Y&F	10	68013
Hygicult™ E/Beta-Gur	10	68267

## Clean Card® Pro

Detects Protein Residues on Surfaces.

- · Easy to use surface hygiene test pad
- Results in 30 seconds
- Long shelf life; stored at room temperature
- · No instrument or incubation required

ORDERING INFORMATION		
Product	# OF TESTS	ITEM #
Clean Card® Pro 25 (includes water spray bottle)	25	133148
Clean Card® Pro 50	50	133149





## **INDUSTRIAL**

#### Easicult®

Convenient Industrial Culture System Identifying Bacteria & Fungi in Oils, Water, Fuels and Other Applications.

- Ready to use on site test for monitoring microbial contamination
- Easy result interpretation
- · Decreases costly disposal of unusable fluids
- Incubate at 36°C for up to 24 hours
- Option to incubate at room temperature for up to 5 days



ORDERING INFORMATION		
Product	# OF TESTS	ITEM#
Easicult® TTC	10	4030
Easicult® Combi	10	4060





## **INCUBATOR**



#### Cultura® M Incubator

Small, Durable and Reliable Incubator. Suitable for the Incubation of Dip Paddles.

- · Lightweight and portable
- Temperature range from 25°C to 48°C
- Holds 18 dip paddles

ORDERING INFORMATION		
Product	# OF TESTS	ITEM #
Cultura® M Incubator	NA	05866

# lifeSign

EARLY DETECTION . RAPID DIAGNOSIS . BETTER OUTCOMES



For pricing and purchasing, please reach out and contact:



**2**48-781-8700



**1-877-745-9943** 



info@elproviders.com



www.elproviders.com



## A PBM Group Company

85 Orchard Road • Skillman • New Jersey • 08558 800.526.2125 • 732.246.3366 • Fax: 732.246.0570