

Automatic Viral DNA/RNA Extraction Using **ExiPrep™ 96Lite**

01

MINIMALISM

Compact size
allows high
space-efficiency

02

EASY

Minimized
hands-on step &
innate protocols

03

SPEEDY

Viral DNA/RNA
extraction of up to 96
samples in 50 minutes



ExiPrep™96Lite




ExiPrep™96 Viral DNA/RNA Kit

ExiPrep™ 96Lite

- **Kit:** ExiPrep™96 Viral DNA/RNA Kit (K-4614)
- **Specimen:** Sputum, Nasopharyngeal & Oropharyngeal swab
- Starting volume : 200-400 $\mu\ell$
- Elution volume : 100 $\mu\ell$
- **Throughput :** 96 samples/run
- **Turnaround Time :** 50min



■ RNA extraction protocol for COVID-19 diagnosis



Ordering Information

Category	Product	Cat. No.	Tests/kit
Extraction Kit	ExiPrep™96 Viral DNA/RNA Kit	K-4614	384 tests
Instrument	ExiPrep™96 Lite	A-5250	-
Diagnostics Kit	AccuPower® SARS-CoV-2 Multiplex Real-Time RT-PCR Kit	SCVM-2112	100 tests
	AccuPower® RV1 Multiplex Kit	RV1-2112	100 tests